

Recombinant Human TRIM37 Protein, N-His

Catalog No: #AP97544

Package Size: #AP97544-1 50ug #AP97544-2 100ug #AP97544-3 1mg

Description

Product Name	Recombinant Human TRIM37 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E3 ubiquitin-protein ligase TRIM37, 2.3.2.27, Mulibrey nanism protein, RING-type E3 ubiquitin transferase TRIM37, Tripartite motif-containing protein 37, TRIM37, KIAA0898, MUL, POB1
Accession No.	O94972
Target Species	Homo sapiens (Human)
Calculated MW	17.44 kDa
Target Length	Ser275-Ser405
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.