

Recombinant Human GCNT3 Protein, N-His

Catalog No: #AP97554

Package Size: #AP97554-1 50ug #AP97554-2 100ug #AP97554-3 1mg

Description

Product Name	Recombinant Human GCNT3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Beta-1,3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 3, 2.4.1.102, 2.4.1.148, 2.4.1.150, C2GnT-mucin type, C2GnT-M, hC2GnT-M, Core 2/core 4 beta-1,6-N-acetylglucosaminyltransferase, C2/4GnT, GCNT3
Accession No.	O95395
Target Species	Homo sapiens (Human)
Calculated MW	50.06 kDa
Target Length	Leu27-Leu438
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.