

Recombinant Human AHSA1 Protein, N-His

Catalog No: #AP97557

Package Size: #AP97557-1 50ug #AP97557-2 100ug #AP97557-3 1mg

Description

Product Name	Recombinant Human AHSA1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Activator of 90 kDa heat shock protein ATPase homolog 1, AHA1, p38, AHSA1, C14orf3, HSPC322
Accession No.	O95433
Target Species	Homo sapiens (Human)
Calculated MW	17.98 kDa
Target Length	His26-Ile164
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.