

Recombinant Human POLR1A Protein, N-His

Catalog No: #AP97562

Package Size: #AP97562-1 50ug #AP97562-2 100ug #AP97562-3 1mg

Description

Product Name	Recombinant Human POLR1A Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	DNA-directed RNA polymerase I subunit RPA1, RNA polymerase I subunit A1, 2.7.7.6, A190, DNA-directed RNA polymerase I largest subunit, DNA-directed RNA polymerase I subunit A, RNA polymerase I 194 kDa subunit, RPA194, POLR1A
Accession No.	O95602
Target Species	Homo sapiens (Human)
Calculated MW	23.40 kDa
Target Length	Asp438-Ala628
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.