

Recombinant Human ESRRB Protein, N-His

Catalog No: #AP97566

Package Size: #AP97566-1 50ug #AP97566-2 100ug #AP97566-3 1mg

Description

Product Name	Recombinant Human ESRRB Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ERR-beta, NR3B2, Estrogen-related receptor beta, ESRL2, Estrogen receptor-like 2, ERR beta-2, ESRRB, Nuclear receptor subfamily 3 group B member 2, Steroid hormone receptor ERR2, ERRB2
Accession No.	O95718
Target Species	Homo sapiens (Human)
Calculated MW	25.21 kDa
Target Length	Gly235-Val433
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.