

Recombinant Human NUDT14 Protein, N-His

Catalog No: #AP97569

Package Size: #AP97569-1 50ug #AP97569-2 100ug #AP97569-3 1mg

Description

Product Name	Recombinant Human NUDT14 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Uridine diphosphate glucose pyrophosphatase NUDT14, UDPG pyrophosphatase, UGPPase, 3.6.1.45, Nucleoside diphosphate-linked moiety X motif 14, Nudix motif 14, NUDT14, UGPP
Accession No.	O95848
Target Species	Homo sapiens (Human)
Calculated MW	26.28 kDa
Target Length	Met1-Gln222
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.