

Recombinant Human BPNT1 Protein, N-His

Catalog No: #AP97572

Package Size: #AP97572-1 50ug #AP97572-2 100ug #AP97572-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BPNT1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	3'(2'),5'-bisphosphate nucleotidase 1, 3.1.3.7, 3'-phosphoadenosine 5'-phosphate phosphatase, PAP phosphatase, Bisphosphate 3'-nucleotidase 1, BPntase 1, HsPIP, Inositol-polyphosphate 1-phosphatase, 3.1.3.57, BPNT1
Accession No.	O95861
Target Species	Homo sapiens (Human)
Calculated MW	35.56 kDa
Target Length	Met1-Pro308
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.