

Recombinant Human TGM6 Protein, N-His

Catalog No: #AP97576

Package Size: #AP97576-1 50ug #AP97576-2 100ug #AP97576-3 1mg

Description

Product Name	Recombinant Human TGM6 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Protein-glutamine gamma-glutamyltransferase 6, 2.3.2.13, Transglutaminase Y, TGY, TGase Y, Transglutaminase-3-like, TGase-3-like, Transglutaminase-6, TG6, TGase-6, TGM6, TGM3L
Accession No.	O95932
Target Species	Homo sapiens (Human)
Calculated MW	25.70 kDa
Target Length	Ser494-Lys706
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.