

Recombinant Human NUDT3 Protein, N-His

Catalog No: #AP97581

Package Size: #AP97581-1 50ug #AP97581-2 100ug #AP97581-3 1mg

Description

Product Name	Recombinant Human NUDT3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Diphosphoinositol polyphosphate phosphohydrolase 1, DIPPP-1, 3.6.1.52, Diadenosine hexaphosphate hydrolase, Ap6A hydrolase, 3.6.1.61, Endopolyphosphatase, 3.6.1.10, Nucleoside diphosphate-linked moiety X motif 3, Nudix motif 3, m7GpppN-mRNA hydrolase, 3.6.1.62, m7GpppX diphosphatase, 3.6.1.59, NUDT3, DIPPP, DIPPP1
Accession No.	O95989
Target Species	Homo sapiens (Human)
Calculated MW	19.08 kDa
Target Length	Gly16-Ser163
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.