

Recombinant Human CYB5A Protein, N-GST & C-His

Catalog No: #AP97583

Package Size: #AP97583-1 50ug #AP97583-2 100ug #AP97583-3 1mg

Description

Product Name	Recombinant Human CYB5A Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Cytochrome b5, Microsomal cytochrome b5 type A, MCB5, CYB5A, CYB5
Accession No.	P00167
Target Species	Homo sapiens (Human)
Calculated MW	37.94 kDa
Target Length	Tyr11-Lys94
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.