

Recombinant Human MOS Protein, N-His

Catalog No: #AP97584

Package Size: #AP97584-1 50ug #AP97584-2 100ug #AP97584-3 1mg

Description

Product Name	Recombinant Human MOS Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Proto-oncogene serine/threonine-protein kinase mos, 2.7.11.1, Oocyte maturation factor mos, Proto-oncogene c-Mos, MOS
Accession No.	P00540
Target Species	Homo sapiens (Human)
Calculated MW	34.29 kDa
Target Length	Ser54-Gly346
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.