

Recombinant Human PPY Protein, N-GST & C-His

Catalog No: #AP97586

Package Size: #AP97586-1 50ug #AP97586-2 100ug #AP97586-3 1mg

Description

Product Name	Recombinant Human PPY Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Pancreatic polypeptide prohormone, PH, Pancreatic polypeptide Y, Obinepitide, Pancreatic polypeptide, HPP, PP, Pancreatic icosapeptide, PI, PPY, PNP
Accession No.	P01298
Target Species	Homo sapiens (Human)
Calculated MW	34.92 kDa
Target Length	Ala30-Arg88
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.