

Recombinant Human PRH1 Protein, N-His-KSI

Catalog No: #AP97590

Package Size: #AP97590-1 50ug #AP97590-2 100ug #AP97590-3 1mg

Description

Product Name	Recombinant Human PRH1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Salivary acidic proline-rich phosphoprotein 1/2, Db-s, PRP-1/PRP-2, Parotid acidic protein, Pa, Parotid double-band protein, Parotid isoelectric focusing variant protein, PIF-S, Parotid proline-rich protein 1/2, Pr1/Pr2, Protein C, Salivary acidic proline-rich phosphoprotein 1/2, Salivary acidic proline-rich phosphoprotein 3/4, Db-F, PIF-F, PRP-3/PRP-4, Protein A, Peptide P-C, PRH1, PRH2
Accession No.	P02810
Target Species	Homo sapiens (Human)
Calculated MW	31.13 kDa
Target Length	Gln17-Gln166
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.