

Recombinant Human MT-ND2 Protein, N-His-KSI

Catalog No: #AP97591

Package Size: #AP97591-1 50ug #AP97591-2 100ug #AP97591-3 1mg

Description

Product Name	Recombinant Human MT-ND2 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NADH-ubiquinone oxidoreductase chain 2, 7.1.1.2, NADH dehydrogenase subunit 2, MT-ND2, MTND2, NADH2, ND2
Accession No.	P03891
Target Species	Homo sapiens (Human)
Calculated MW	25.81 kDa
Target Length	[Trp30-Lys58] & [Pro116-Ser148] & [Ser299-Thr323]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.