

Recombinant Human F13B Protein, N-His

Catalog No: #AP97593

Package Size: #AP97593-1 50ug #AP97593-2 100ug #AP97593-3 1mg

Description

Product Name	Recombinant Human F13B Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Coagulation factor XIII B chain, Fibrin-stabilizing factor B subunit, Protein-glutamine gamma-glutamyltransferase B chain, Transglutaminase B chain, F13B
Accession No.	P05160
Target Species	Homo sapiens (Human)
Calculated MW	38.42 kDa
Target Length	Lys333-Gln650
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.