

Recombinant Human RPLP0 Protein, N-His

Catalog No: #AP97596

Package Size: #AP97596-1 50ug #AP97596-2 100ug #AP97596-3 1mg

Description

Product Name	Recombinant Human RPLP0 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	60S acidic ribosomal protein P0, 60S ribosomal protein L10E, Large ribosomal subunit protein uL10, RPLP0
Accession No.	P05388
Target Species	Homo sapiens (Human)
Calculated MW	35.37 kDa
Target Length	Pro2-Asp317
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.