

Recombinant Human POLR3D Protein, N-His

Catalog No: #AP97597

Package Size: #AP97597-1 50ug #AP97597-2 100ug #AP97597-3 1mg

Description

Product Name	Recombinant Human POLR3D Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	DNA-directed RNA polymerase III subunit RPC4, RNA polymerase III subunit C4, DNA-directed RNA polymerase III subunit D, Protein BN51, RNA polymerase III 47 kDa subunit, RPC53 homolog, POLR3D, BN51, BN51T
Accession No.	P05423
Target Species	Homo sapiens (Human)
Calculated MW	18.75 kDa
Target Length	Val249-Arg398
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.