

Recombinant Human ADH4 Protein, N-His

Catalog No: #AP97609

Package Size: #AP97609-1 50ug #AP97609-2 100ug #AP97609-3 1mg

Description

Product Name	Recombinant Human ADH4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	All-trans-retinol dehydrogenase [NAD(+)] ADH4, 1.1.1.105, Alcohol dehydrogenase 2, Alcohol dehydrogenase 4, Alcohol dehydrogenase class II pi chain, ADH4, ADH2
Accession No.	P08319
Target Species	Homo sapiens (Human)
Calculated MW	42.38 kDa
Target Length	Met1-Phe380
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.