

Recombinant Human SNRPA Protein, N-His

Catalog No: #AP97611

Package Size: #AP97611-1 50ug #AP97611-2 100ug #AP97611-3 1mg

Description

Product Name	Recombinant Human SNRPA Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	U1 snRNP A, U1 small nuclear ribonucleoprotein A, SNRPA, U1-A, U1A
Accession No.	P09012
Target Species	Homo sapiens (Human)
Calculated MW	13.89 kDa
Target Length	Gly99-His208
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.