

Recombinant Human TFPT Protein, N-GST & C-His

Catalog No: #AP97617

Package Size: #AP97617-1 50ug #AP97617-2 100ug #AP97617-3 1mg

Description

Product Name	Recombinant Human TFPT Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	TCF3 fusion partner, INO80 complex subunit F, Protein FB1, TFPT, INO80F
Accession No.	P0C1Z6
Target Species	Homo sapiens (Human)
Calculated MW	37.57 kDa
Target Length	Gly53-Gly128
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.