

Recombinant Human POLR1D Protein, N-His

Catalog No: #AP97622

Package Size: #AP97622-1 50ug #AP97622-2 100ug #AP97622-3 1mg

Description

Product Name	Recombinant Human POLR1D Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	DNA-directed RNA polymerases I and III subunit RPAC2, RNA polymerases I and III subunit AC2, AC19, DNA-directed RNA polymerase I subunit D, RNA polymerase I 16 kDa subunit, RPA16, RPC16, hRPA19, POLR1D
Accession No.	P0DPB6
Target Species	Homo sapiens (Human)
Calculated MW	17.40 kDa
Target Length	Met1-Phe133
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.