

Recombinant Human UROS Protein, N-His

Catalog No: #AP97626

Package Size: #AP97626-1 50ug #AP97626-2 100ug #AP97626-3 1mg

Description

Product Name	Recombinant Human UROS Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Uroporphyrinogen-III synthase, UROIIS, UROS, 4.2.1.75, Hydroxymethylbilane hydrolyase [cyclizing], Uroporphyrinogen-III cosynthase, UROS
Accession No.	P10746
Target Species	Homo sapiens (Human)
Calculated MW	30.79 kDa
Target Length	Met1-Cys265
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.