

Recombinant Human SLC2A1 Protein, N-His-KSI

Catalog No: #AP97628

Package Size: #AP97628-1 50ug #AP97628-2 100ug #AP97628-3 1mg

Description

Product Name	Recombinant Human SLC2A1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Solute carrier family 2, facilitated glucose transporter member 1, Glucose transporter type 1, erythrocyte/brain, GLUT-1, HepG2 glucose transporter, SLC2A1, GLUT1
Accession No.	P11166
Target Species	Homo sapiens (Human)
Calculated MW	32.20 kDa
Target Length	[Asn34-Ser66] & [Cys207-Pro271] & [Lys451-Val492]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.