

Recombinant Human NAT2 Protein, N-His

Catalog No: #AP97631

Package Size: #AP97631-1 50ug #AP97631-2 100ug #AP97631-3 1mg

Description

Product Name	Recombinant Human NAT2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Arylamine N-acetyltransferase 2, 2.3.1.5, Arylamide acetylase 2, N-acetyltransferase type 2, NAT-2, Polymorphic arylamine N-acetyltransferase, PNAT, NAT2, AAC2
Accession No.	P11245
Target Species	Homo sapiens (Human)
Calculated MW	35.73 kDa
Target Length	Met1-Ile290
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.