

Recombinant Human RPS2 Protein, N-His

Catalog No: #AP97645

Package Size: #AP97645-1 50ug #AP97645-2 100ug #AP97645-3 1mg

Description

Product Name	Recombinant Human RPS2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Small ribosomal subunit protein uS5, 40S ribosomal protein S2, 40S ribosomal protein S4, Protein LLRep3, RPS2, RPS4
Accession No.	P15880
Target Species	Homo sapiens (Human)
Calculated MW	22.60 kDa
Target Length	Thr64-Tyr250
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.