

Recombinant Human SLC6A1 Protein, N-His-KSI

Catalog No: #AP97667

Package Size: #AP97667-1 50ug #AP97667-2 100ug #AP97667-3 1mg

Description

Product Name	Recombinant Human SLC6A1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Sodium- and chloride-dependent GABA transporter 1, GAT-1, Solute carrier family 6 member 1, SLC6A1, GABATR, GABT1, GAT1
Accession No.	P30531
Target Species	Homo sapiens (Human)
Calculated MW	23.77 kDa
Target Length	Ser145-Arg211
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.