

Recombinant Human DLST/OGDC-E2 Protein, N-His-KSI

Catalog No: #AP97674

Package Size: #AP97674-1 50ug #AP97674-2 100ug #AP97674-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DLST/OGDC-E2 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Dihydrolipoyllysine-residue succinyltransferase component of 2-oxoglutarate dehydrogenase complex, mitochondrial, 2.3.1.61, 2-oxoglutarate dehydrogenase complex component E2, OGDC-E2, Dihydroipoamide succinyltransferase component of 2-oxoglutarate dehydrogenase complex, E2K, DLST, DLTS
Accession No.	P36957
Target Species	Homo sapiens (Human)
Calculated MW	39.02 kDa
Target Length	Leu2-Lys226
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.