

Recombinant Human CCT4 Protein, N-His

Catalog No: #AP97681

Package Size: #AP97681-1 50ug #AP97681-2 100ug #AP97681-3 1mg

Description

Product Name	Recombinant Human CCT4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	CCTD, T-complex protein 1 subunit delta, TCP-1-delta, Stimulator of TAR RNA-binding, CCT4, CCT-delta, SRB
Accession No.	P50991
Target Species	Homo sapiens (Human)
Calculated MW	57.99 kDa
Target Length	Gln25-Arg539
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.