

## Recombinant Human INSL3 Protein, N-GST &amp; C-His

Catalog No: #AP97682

Package Size: #AP97682-1 50ug #AP97682-2 100ug #AP97682-3 1mg

## Description

Product Name	Recombinant Human INSL3 Protein, N-GST & C-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Insulin-like 3, Leydig insulin-like peptide, Ley-I-L, Relaxin-like factor, Insulin-like 3 B chain, Insulin-like 3 A chain, INSL3, RLF, RLNL
Accession No.	P51460
Target Species	Homo sapiens (Human)
Calculated MW	40.55 kDa
Target Length	Leu21-Tyr131
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.