

Recombinant Human DEFB1 Protein, N-His-KSI

Catalog No: #AP97689

Package Size: #AP97689-1 50ug #AP97689-2 100ug #AP97689-3 1mg

Description

Product Name	Recombinant Human DEFB1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Beta-defensin 1, BD-1, hBD-1, Defensin, beta 1, DEFB1, BD1, HBD1
Accession No.	P60022
Target Species	Homo sapiens (Human)
Calculated MW	20.84 kDa
Target Length	Gly22-Lys68
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.