

Recombinant Human SLC6A3 Protein, N-His-KSI

Catalog No: #AP97693

Package Size: #AP97693-1 50ug #AP97693-2 100ug #AP97693-3 1mg

Description

Product Name	Recombinant Human SLC6A3 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Sodium-dependent dopamine transporter, DA transporter, DAT, Solute carrier family 6 member 3, SLC6A3, DAT1
Accession No.	Q01959
Target Species	Homo sapiens (Human)
Calculated MW	31.93 kDa
Target Length	[His165-Pro235] & [Thr286-Ser309] & [Gln373-Ala399] & [Leu474-His477] & [Thr542-Ala557]
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.