

Recombinant Human GUCY1A1 Protein, N-His

Catalog No: #AP97694

Package Size: #AP97694-1 50ug #AP97694-2 100ug #AP97694-3 1mg

Description

Product Name	Recombinant Human GUCY1A1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Guanylate cyclase soluble subunit alpha-1, GCS-alpha-1, 4.6.1.2, Guanylate cyclase soluble subunit alpha-3, GCS-alpha-3, Soluble guanylate cyclase large subunit, GUCY1A1, GUC1A3, GU CSA3, GUCY1A3
Accession No.	Q02108
Target Species	Homo sapiens (Human)
Calculated MW	26.76 kDa
Target Length	Gln468-Asp690
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.