

Recombinant Human USP4 Protein, N-His

Catalog No: #AP97698

Package Size: #AP97698-1 50ug #AP97698-2 100ug #AP97698-3 1mg

Description

Product Name	Recombinant Human USP4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ubiquitin carboxyl-terminal hydrolase 4, 3.4.19.12, Deubiquitinating enzyme 4, Ubiquitin thioesterase 4, Ubiquitin-specific-processing protease 4, Ubiquitous nuclear protein homolog, USP4, UNP, UNPH
Accession No.	Q13107
Target Species	Homo sapiens (Human)
Calculated MW	26.67 kDa
Target Length	Asp12-Glu220
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.