

Recombinant Human NEUROD1 Protein, N-His-KSI

Catalog No: #AP97699

Package Size: #AP97699-1 50ug #AP97699-2 100ug #AP97699-3 1mg

Description

Product Name	Recombinant Human NEUROD1 Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Neurogenic differentiation factor 1, NeuroD, NeuroD1, Class A basic helix-loop-helix protein 3, bHLHa3, NEUROD1, BHLHA3, NEUROD
Accession No.	Q13562
Target Species	Homo sapiens (Human)
Calculated MW	55.55 kDa
Target Length	Thr2-Asp356
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.