

Recombinant Human TARBP2 Protein, N-His

Catalog No: #AP97705

Package Size: #AP97705-1 50ug #AP97705-2 100ug #AP97705-3 1mg

Description

Product Name	Recombinant Human TARBP2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	RISC-loading complex subunit TARBP2, TAR RNA-binding protein 2, Trans-activation-responsive RNA-binding protein, TARBP2, TRBP
Accession No.	Q15633
Target Species	Homo sapiens (Human)
Calculated MW	13.30 kDa
Target Length	Thr264-Lys366
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.