

Recombinant Human INF2 Protein, N-His

Catalog No: #AP97707

Package Size: #AP97707-1 50ug #AP97707-2 100ug #AP97707-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human INF2 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Inverted formin-2, HBEBP2-binding protein C, INF2, C14orf151, C14orf173
Accession No.	Q27J81
Target Species	Homo sapiens (Human)
Calculated MW	35.56 kDa
Target Length	Glu654-Glu946
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.