

Recombinant Human ACSM3 Protein, N-His

Catalog No: #AP97709

Package Size: #AP97709-1 50ug #AP97709-2 100ug #AP97709-3 1mg

Description

Product Name	Recombinant Human ACSM3 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Acyl-coenzyme A synthetase ACSM3, mitochondrial, 6.2.1.2, Acyl-CoA synthetase medium-chain family member 3, Butyrate--CoA ligase 3, Butyryl-coenzyme A synthetase 3, Middle-chain acyl-CoA synthetase 3, Propionate--CoA ligase, 6.2.1.17, Protein SA homolog, ACSM3, SAH
Accession No.	Q53FZ2
Target Species	Homo sapiens (Human)
Calculated MW	33.22 kDa
Target Length	His302-Glu577
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.