

Recombinant Human HSD17B12 Protein, N-His

Catalog No: #AP97710

Package Size: #AP97710-1 50ug #AP97710-2 100ug #AP97710-3 1mg

Description

Product Name	Recombinant Human HSD17B12 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Very-long-chain 3-oxoacyl-CoA reductase, 1.1.1.330, 17-beta-hydroxysteroid dehydrogenase 12, 17-beta-HSD 12, 3-ketoacyl-CoA reductase, KAR, Estradiol 17-beta-dehydrogenase 12, 1.1.1.62, Short chain dehydrogenase/reductase family 12C member 1, HSD17B12, SDR12C1
Accession No.	Q53GQ0
Target Species	Homo sapiens (Human)
Calculated MW	19.29 kDa
Target Length	Ser27-Lys181
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.