

Recombinant Human CCDC80 Protein, N-His

Catalog No: #AP97711

Package Size: #AP97711-1 50ug #AP97711-2 100ug #AP97711-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CCDC80 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Coiled-coil domain-containing protein 80, Down-regulated by oncogenes protein 1, Up-regulated in BRS-3 deficient mouse homolog, CCDC80, DRO1, URB, HBE245
Accession No.	Q76M96
Target Species	Homo sapiens (Human)
Calculated MW	20.21 kDa
Target Length	Asn140-Gly297
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.