

Recombinant Human PDPN/Podoplanin Protein, N-His-KSI

Catalog No: #AP97712

Package Size: #AP97712-1 50ug #AP97712-2 100ug #AP97712-3 1mg

Description

Product Name	Recombinant Human PDPN/Podoplanin Protein, N-His-KSI
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Podoplanin, Aggrus, Glycoprotein 36, Gp36, PA2.26 antigen, T1-alpha, T1A, 29kDa cytosolic podoplanin intracellular domain, PICD, PDPN, GP36, PSEC0003, PSEC0025
Accession No.	Q86YL7
Target Species	Homo sapiens (Human)
Calculated MW	26.89 kDa
Target Length	Ala23-Leu131
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.