

Recombinant Human NOX1 Protein, N-His

Catalog No: #AP97717

Package Size: #AP97717-1 50ug #AP97717-2 100ug #AP97717-3 1mg

Description

Product Name	Recombinant Human NOX1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	NADPH oxidase 1, NOX-1, 1.6.3.-, Mitogenic oxidase 1, MOX-1, NADH/NADPH mitogenic oxidase subunit P65-MOX, NOH-1, NOX1, MOX1, NOH1
Accession No.	Q9Y5S8
Target Species	Homo sapiens (Human)
Calculated MW	19.48 kDa
Target Length	Trp418-Phe564
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.