

Recombinant Human RNF168 Protein, N-His

Catalog No: #AP97719

Package Size: #AP97719-1 50ug #AP97719-2 100ug #AP97719-3 1mg

Description

Product Name	Recombinant Human RNF168 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	E3 ubiquitin-protein ligase RNF168, hRNF168, 2.3.2.27, RING finger protein 168, RING-type E3 ubiquitin transferase RNF168, RNF168
Accession No.	Q8IYW5
Target Species	Homo sapiens (Human)
Calculated MW	13.08 kDa
Target Length	Met1-Gln95
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.