

Recombinant Human ERVW-1/HERV-W-Env Protein, N-His

Catalog No: #AP97720

Package Size: #AP97720-1 50ug #AP97720-2 100ug #AP97720-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ERVW-1/HERV-W-Env Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Syncytin-1, Endogenous retrovirus group W member 1, Syncytin, env
Accession No.	Q96TB1
Target Species	Homo sapiens (Human)
Calculated MW	22.05 kDa
Target Length	Asn1-Leu180
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.