

# Recombinant Human SCARB1/SR-BI/CLA-1/SRB1 Protein, N-His

Catalog No: #AP97721

Package Size: #AP97721-1 50ug #AP97721-2 100ug #AP97721-3 1mg

## Description

|                |                                                                                                                                                                                                 |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name   | Recombinant Human SCARB1/SR-BI/CLA-1/SRB1 Protein, N-His                                                                                                                                        |
| Host Species   | E. coli                                                                                                                                                                                         |
| Purification   | >90% as determined by SDS-PAGE.                                                                                                                                                                 |
| Applications   | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress                                                                                                                                 |
| Other Names    | Scavenger receptor class B member 1, SRB1, CD36 and LIMPII analogous 1, CLA-1, CD36 antigen-like 1, Collagen type I receptor, thrombospondin receptor-like 1, SR-BI, CD36, SCARB1, CD36L1, CLA1 |
| Accession No.  | Q8WTV0                                                                                                                                                                                          |
| Target Species | Homo sapiens (Human)                                                                                                                                                                            |
| Calculated MW  | 43.38 kDa                                                                                                                                                                                       |
| Target Length  | Ile61-Gly424                                                                                                                                                                                    |
| Tag Info       | His-tag                                                                                                                                                                                         |
| Formulation    | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.                                                                                                      |
| Storage        | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.                       |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.