

Recombinant Human PTP4A1/PRL1 Protein, N-His

Catalog No: #AP97723

Package Size: #AP97723-1 50ug #AP97723-2 100ug #AP97723-3 1mg

Description

Product Name	Recombinant Human PTP4A1/PRL1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	PRL-1, PTP(CAAXI), PTPCAAX1, PTP4A1, Protein tyrosine phosphatase type IVA 1, PRL1, Protein-tyrosine phosphatase of regenerating liver 1, Protein-tyrosine phosphatase 4a1
Accession No.	Q93096
Target Species	Homo sapiens (Human)
Calculated MW	21.98 kDa
Target Length	Met1-Gln173
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.