

Recombinant Human IL4I1 Protein, N-His

Catalog No: #AP97727

Package Size: #AP97727-1 50ug #AP97727-2 100ug #AP97727-3 1mg

Description

Product Name	Recombinant Human IL4I1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	L-amino-acid oxidase, LAAO, LAO, 1.4.3.2, 1.4.3.25, Interleukin-4-induced protein 1, IL4-induced protein 1, hIL4I1, Protein Fig-1, hFIG1, IL4I1, FIG1, UNQ636/PRO1265
Accession No.	Q96RQ9
Target Species	Homo sapiens (Human)
Calculated MW	56.76 kDa
Target Length	Gln22-Asn505
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.