

Recombinant Human LRRC4 Protein, N-His

Catalog No: #AP97732

Package Size: #AP97732-1 50ug #AP97732-2 100ug #AP97732-3 1mg

Description

Product Name	Recombinant Human LRRC4 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Leucine-rich repeat-containing protein 4, Brain tumor-associated protein BAG, Nasopharyngeal carcinoma-associated gene 14 protein, Netrin-G2 ligand, NGL-2, LRRC4, BAG, NAG14, UNQ554/PRO1111
Accession No.	Q9HBW1
Target Species	Homo sapiens (Human)
Calculated MW	47.82 kDa
Target Length	Ala39-Thr443
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.