

Recombinant Human USP31 Protein, N-His

Catalog No: #AP97737

Package Size: #AP97737-1 50ug #AP97737-2 100ug #AP97737-3 1mg

Description

Product Name	Recombinant Human USP31 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Ubiquitin carboxyl-terminal hydrolase 31, 3.4.19.12, Deubiquitinating enzyme 31, Ubiquitin thioesterase 31, Ubiquitin-specific-processing protease 31, USP31, KIAA1203
Accession No.	Q70CQ4
Target Species	Homo sapiens (Human)
Calculated MW	27.56 kDa
Target Length	Thr595-Ser816
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.