

Recombinant Human SUPV3L1 Protein, N-His

Catalog No: #AP97738

Package Size: #AP97738-1 50ug #AP97738-2 100ug #AP97738-3 1mg

Description

Product Name	Recombinant Human SUPV3L1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	ATP-dependent RNA helicase SUPV3L1, mitochondrial, 3.6.4.13, Suppressor of var1 3-like protein 1, SUV3-like protein 1, SUPV3L1, SUV3
Accession No.	Q8IYB8
Target Species	Homo sapiens (Human)
Calculated MW	78.21 kDa
Target Length	Thr47-Glu722
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.