

Recombinant Human SUN5 Protein, N-His

Catalog No: #AP97743

Package Size: #AP97743-1 50ug #AP97743-2 100ug #AP97743-3 1mg

Description

Product Name	Recombinant Human SUN5 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	SUN domain-containing protein 5, Sad1 and UNC84 domain-containing protein 5, Sperm-associated antigen 4-like protein, Testis and spermatogenesis-related gene 4 protein, SUN5, SPAG4L, TSARG4
Accession No.	Q8TC36
Target Species	Homo sapiens (Human)
Calculated MW	20.44 kDa
Target Length	Lys202-Ser362
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.